

Endoscopy Cleaning Brush

Single-Use Brush

- A single-use brush for use in association with water and other solvents to remove organic soils and other residuals from channels and valves of flexible endoscopes
- The endoscopy brush is sufficiently flexible for the smallest of channels
- Brushes are single-use only and cannot be re-used or re-processed
- Comply with the essential requirements of the Medical Devices Directive EC93/42-2007/47/EC - Class 1
- They are CE marked on the product label



Technical Details

Material	Nylon bristles and tubing
Dimensions	Bristle Length: 9 mm
	Maximum Diameter: 7.5 mm
	Tubing Length: 2300 mm
	Tubing Outside Diameter: 1.8 mm

Order Code	Description
CBBR005	Clinipak Endoscopy Cleaning Brush