

Flexistem Microclene

Flexible cleaning of small bore items

A development of Flexistem, Microclene is more brush like and finer. FLEX021 will clean down to very small diameters and is ideal for Zoellner suction tubes.

Microclene is a fine bore cleaning system for instruments with a very small bore. It can be used to safely dislodge deposits during the washing process.

Flexistem Microclene is used as a brush ensuring that the final rinse process is very effective. As with all Flexistem products, Microclene can be doubled over to increase the cleaning diameter or to increase the rigidity of the wire. Supplied in lengths of 7.5 metre rolls which can be cut to length with the cutters provided, Microclene may be safely disposed of after use.



Order Code	Description
CFLEX021	Clinipak Flexistem Microclene 3mm Super Fine Brush
CFLEX022	Clinipak Flexistem Microclene 6mm Medium Fine Brush
CFLEX023	Clinipak Flexistem Microclene 3mm Super Brush
CFLEX024	Clinipak Flexistem Microclene 6mm Medium Brush
CFLEX025	Clinipak Flexistem Microclene Complete set including dispenser, cutter & one pack of each item
CFLEX026	Clinipak Flexistem Microclene Refill Set (one pack of each item)

